fms solo[®]


Fluid Management System

Instruction Manual

CE 0197

fms solo® system

English

A. 2 portals in the knee

Español

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Italiano

A. 2 accessi al ginocchio

Deutsch

A. 2 Zugänge im Knie

Français

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Português

A. 2 portas no joelho

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Unpacking, assembly and initial system check

Unpacking:

Verify that all items on packing list match the products you have received (verify serial numbers and references). Damage should be reported immediately to FMS Customer Service

Department: Tel: +33 4 92 12 33 93 or
email: repair@fmsgroup.com

Save all packing material, as they will be required to ship products in the event of a product return.

Assembly:

- Connect the power cord to the power connection (14).
- Connect the foot pedal to the foot pedal connection (16) on the back of the pump.

Initial System Check:

- Turn the device ON by pressing the POWER (4) pad.
- Switch the I/O switch (14) on "I".
- Press the Pump Run/Stop pads (10) to Run.
- Press on all pads (an audible tone should be heard for each pad pressed).

Note: numbers in (), refer to diagram in section 1.6.

Weekly inspection:

FMS recommends a weekly inspection of all equipment. These inspections can be performed on a less regular basis depending on the duration and frequency of use.

- Perform a visual inspection of the device. If damage is noted return to FMS Service Centre for internal damage inspection.
- Inspect that the transparent safety cover (1) and black hinges are in good working order.
- Inspect power cord and all cables for wear or cuts.
- Ensure that the foot pedal and shaver interface connector pins are not missing or bent.
- Plug the power cord to the power connection.
- Switch the I/O switch (14) on "I".
- Turn the device ON by pressing the **POWER** (4) pad.
- Press the Pump Run/Stop pads (10).
- Press on all pads (an audible tone should be heard for each pad pressed).

Annual inspection:

FMS recommends that a bio-medical engineer or authorized service personal annually check all equipment to assess its functionality and technical safety.

Refer to Service manual for annual inspection procedures.

Table of contents

CHAPTER 1. INTRODUCTION	8
1.1. Product description	8
1.2. Declaration of conformity	8
1.3. General warnings and cautions	8
1.4. Intended use	9
1.5. Symbol definitions	10
1.6. System indicators	12
A. Front panel	12
B. Rear panel	14
1.7. Other indicators	14
CHAPTER 2. OPERATING PRECAUTIONS, CHECKS & WARNINGS	15
2.1. Sterile packaging	15
2.2. Storage and handling	15
CHAPTER 3. SET-UP PROCEDURES	16
3.1. Pump set-up	16
3.2. Tubing set-up	17
A. Installing the “One Day Set™” Ref. 4503	18
B. Irrigation tubing set up Ref. 1102	20
C. Getting Started	21
D. End of operation	22
E. Set up of next operation	23
F. End of operating day	24

CHAPTER 4. PUMP SYSTEM FUNCTIONS	25
4.1. Run/stop	25
4.2. Preset pressure	25
4.3. Filling the chamber	25
4.4. Recommended initial preset pressure	26
4.5. Pressure alarm	26
4.6. Pressure Cycle mode (red pedal)	27
4.7. Maintenance: cleaning & decontamination	27
4.8. Interface cables	28
CHAPTER 5. TROUBLESHOOTING	29
CHAPTER 6. PRODUCT SPECIFICATIONS & ORDERING INFORMATION	31

Chapter 1. Introduction

1.1. Product description

Irrigation system.

The fms solo® is an advanced, and user friendly irrigation system that uses patented technology to improve visualization and lower cost.

Flow compensation.

The fms solo® communicates with your current shaver, improving visualization by activating the “Flow compensator” cycle while the shaver is in rotation.

Pressure cycle red pedal.

One of the most critical aspects of arthroscopy is visualization. The field of view is often obstructed by blood and debris. To alleviate this problem the fms solo® has a pressure cycle. By stepping on the red foot pedal you will activate a timed cycle of increased pressure to control bleeding. Once the cycle is completed, the fms solo® will return to its original configuration.

1.2. Declaration of conformity

We, at fms group declare under our sole responsibility that the product fms solo® is in conformity according to CE 93/42 EEC Medical Device Directive. This product is manufactured in an ISO 9001:2000/13485 facility.

1.3. General warnings and cautions

It is recommended that hospital personnel read this manual before operating, cleaning or sterilizing this product and its accessories. Failure to follow these instructions may result in potential injury and damage or malfunction of the equipment.

The manufacturer and licensed seller of this device do not accept any liability for direct or consequential damage or injury caused by improper use of disposables other than FMS products. Any alterations to this device, repair from an unlicensed service center, or use of non-FMS disposables, will void CE marking, FMS warranty and product liability coverage.

**Warning:**

The fms solo[®] is a Class I Type B device. Ensure that you are using the correct interface cable when connecting to another shaver. For more details see section 4.8. Interface cables.

The system may be affected by electromagnetic interferences from other instruments. Verify that all other instruments and appliances associated or not with the system comply with the standard EN 60601-1-2 (EMC).

If the system continues to be affected, isolate it from the offending instrument and connect to a different main output.

The fms solo[®] is covered by US patent N° 5 131 823, European patent N° 0 448 909 B1, Japanese patent N° 2 892 852. fms solo[®] is a registered trademark of Future Medical System S.A.

Reproduction, transfer, distribution of part or all of the contents in this document in any form without the prior written permission of FMS is prohibited.

1.4. Intended use

The fms solo[®] is intended to be used in arthroscopic surgery in the following joints: shoulder, knee, ankle, elbow, wrist, hip.

1.5. Symbol definitions

Front panel

Power pad

Turns the power ON & OFF if the main power switch on the back of the pump is on (the green standby light is on).

Alarm

The red indicator will light up if a safety parameter is violated.

Pressure cycle.

Alleviates bleeding during surgery, by increasing base pressure for two minutes.

1 – way foot pedal board

Used to activate or cancel a pressure cycle.

Increase

Decrease

Back panel

On

Off

Equipotentiality symbol

IP54

Waterproof for projections

Type B Class 1

Conforms to standard EN 60601-1.

CE mark

and identification of notified body.

Warning

see instructions for use.

Explosion hazard

Main fuse

Replace fuses as marked

Per Directive 2002/96/EC, product must be collected separately. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information.

Disposable Product labeling

Expiration date

Do not wet

Do not reuse

ONE DAY

To be used for only one day

**Date of manufacture
year and month**

**Store between
these temperatures**

Latex free


Sterile

unless the package
is damaged or open.

Sterilization by ethylene oxide.


1.6. System indicators

- 1 Transparent safety covers.**
This cover **must** be closed for the roller pump to rotate.
- 2 IRRIGATION roller pump.**
Provides fluid inflow.
- 3 Tension rocker arm.**
Holds the tubing on the roller pump head.
- 4 POWER pad.**
Turns the power **ON & OFF** if the main power switch in the back of the pump is on (the green “standby” light is on).
- 5 PRESSURE display.**
Displays both the preset pressure when the base pressure is being set, and the dynamic pressure (1 sec) following the release of the pressure adjustment pads.
- 6 ALARM.**
The red indicator will light up if a safety parameter is violated.
- 7 Remote Control connector.**
Connects the remote control to the fms solo®.

- 8 Pressure cycle.**
Alleviates bleeding by increasing base pressure for a 2 minute cycle.
- 9 PRESSURE adjustment pads.**
Press the arrows to increase or decrease the pressure in increments of 5. The initial setting is 50.
- 10 Pump Run/Stop pad.**
Turns the pump ON and OFF. If the pump is activated, <PUMP> appears on the pressure display. Green LED blinks when pump is on stop.
- 11 Built-in Pressure transducer connection.**
To connect the "One Day Set™" and measures the pressure in the pressure chamber.
- 12 Fill Chamber pad.**
Manually fills the chamber and turns off the LOW PRESSURE alarm.
- 13 Patented auto-locking device.**
Allows the disposable tubing to be properly positioned and locked in place around the roller pump.

B – Rear panel

14 I/O switch, fuses and power connection.

15 Identification & specification.

16 Foot pedal connection.

9 - pin connector to plug foot pedal.

17 Shaver interface connection.

25 - pin hand and foot control interface cables connector.

18 Serial number.

- . The first two digits represent the year of manufacturing.

- . The second set of digits represents the week of manufacturing and serial number.

19 Reference number.

1.7. Other indicators (Audible tones).

- An audible tone sounds when the **POWER** pad is pressed and the pump turns on.
- An audible tone sounds when the **PRESSURE** mode is activated (see pressure cycle section to activate the cycle).
- Two audible tones sound when the **PRESSURE** mode is deactivated (see chapter 4 to deactivate the cycle).

Chapter 2. Operating precautions, checks and warnings

Warning:

Electrical safety testing should be performed by a biomedical engineer or other qualified person.

Risk of electrical shock: do not remove the cover. Refer servicing to qualified or authorized fms Repair Service Center.

Do not connect the device to a power source that is not properly earthed (grounded).

Disconnect the device from the main power source when cleaning, servicing, or inspecting.

Inspect all equipment and cables periodically for wear. Replace and return to fms Repair Service Center if damage is noted.

Avoid fluid contact with the fms solo[®], and its electrical connectors.

Do not use flammable agents when cleaning and disinfecting the fms solo[®].

To avoid risk of fire, replace fuses with same type and rating.

2.1. Sterile Packaging:

Carefully examine the shipping package and sterile wrap. If the package is damaged, the sterile seal is broken, or the Expiration Date is passed, do not use. The "One Day Set™" Ref. 4503, can be reused for an entire surgical day not exceeding eight consecutive hours.

Biological and viral studies have shown the efficacy against contamination from one patient to the next during the operating day. Sterility is guaranteed only if the tube changing procedures have been performed precisely according to our set-up guide in Chapter 3.

Studies can be downloaded on our website at:

In English:

www.fmsgroup.com/files/sales/viralen.pdf
www.fmsgroup.com/files/sales/biologicalen.pdf

In French:

www.fmsgroup.com/files/sales/viralfr.pdf
www.fmsgroup.com/files/sales/biologicfr.pdf

All other tube sets, are marked with a  indicating that they are single-use only

2.2. Storage and Handling:

Storage temperature: Do not expose the fms solo[®] to temperatures above 50°C (122°F) or below -10°C (14°F).

Operating Conditions: Do not operate the pump in conditions +10°C (50°F) to +40°C (104°F).

Damage: Do not use a fms solo[®] that has been dropped or shows signs of damage.

Disposal: Do not incinerate. Return to fms for safe disposal.


Factory Settings: Factory settings are based on experience, they should only be used as guidelines. The surgeon is responsible for settings pertaining to his surgical procedure.

Chapter 3. Set-up Procedures

3.1. Pump set-up.

Power on back of the pump on "I".



Open the cover.

Press **Power** pad.
The LCD will light up.



3.2. Tubing set-up.


A. Installing the “One Day Set™” Ref. 4503.

This tubing connects the saline solution bags to the irrigation tubing. Solution goes from the spikes


(A) around the roller pump, through the pressure chamber, where the pressure is read through a filter (B) by a built-in pressure transducer.


Close all three clamps and connect to fluid bags.


Slide the green part of the tubing into the auto-locking device on the pump.


Place and **center** the tubing around the roller pump.

Pull on the tubing and close the tension rocker arm against the tube.

Position the pressure chamber in the holder. **Connect** the transducer to the pump by firmly screwing it in clockwise.

Close the cover.

Green LED blinks. Press **Run/ Stop**. The green LED blinks initially, then remains on.


Press the **Fill chamber** pad several times until the chamber is 1/3 full. The green LED blinks until pad is pressed.


Warning:

It is mandatory that the pressure chamber remains vertical during the procedure and that the pump be positioned at the same level as the patient.


B. Installing the Irrigation tubing Ref. 1102


This tubing connects the irrigation tubing to the sheath in a sterile manner.


After checking that the luer-lock is properly tightened, the sterile nurse takes the irrigation tubing and passes the blue end (D) to a non-sterile person.

The non sterile person removes the blue caps from the irrigation tubing (D) and the "One Day Set™" (C)...


... and **immediately** connects the "One Day Set™" (C) luer-lock, to the irrigation tubing (D) luer-lock.

The pump must be in the **Run** mode (10), with the light on.
Prime the irrigation tubing by opening the clamp under the pressure chamber. When primed, close the clamp or the stop-cock on the arthroscope (in the sterile zone).


C. Getting started.


Insert the cannula tubing onto the cannula.

Insert the shaver tubing onto the shaver suction port.

Connect the luer-lock to the arthroscope. Open the stop-cock on the scope.

Open the clamp under the chamber.


D. End of the operation.

Close the clamp under the pressure chamber.

Press Run/Stop pad (or open transparent safety cover).

Disconnect the blue luer-lock (E) of the irrigation tubing.

Warning:
The sterile buffer zone must remain attached to the "One Day Set™".


Warning: When another surgery is scheduled for the same operating day, leave the "One Day Set™" in place.


The sterile buffer zone must be left on the "One Day Set™", which ensure sterility between surgeries.

E. Set up for next operation.


The "One Day Set™" should have been left in placed at the end of the previous operation with the sterile buffer zone attached.


After checking that the luer-lock is tightened, the sterile nurse takes the irrigation tubing and passes the blue end (D) to a non sterile person.

The non-sterile person removes the old sterile buffer zone by disconnecting the luer-locks of the "One Day Set™" (C) and the irrigation tubing..


... and **immediately** connects the blue luer-lock (D) of the new irrigation tubing to the blue luer-lock (C) of the "One Day Set™".

Warning:
Never attach one
buffer zone
to another.

F. End of operating day.

Close all white clamps.


Turn off the pump by pressing the **POWER** pad. Throw away **all of the tubing.**

Chapter 4. Pump System Functions

4.1. Run/Stop

Press the Run/Stop pad to start all pump system functions.

4.2. Preset pressure

The pressure is displayed in brackets (eg. <50>) on the LCD (5) for 2 seconds when one of the pressure pads is pressed, then the actual pressure is displayed (pressure settings are in levels).

Level 50, equals 5 feet of gravity pressure or the equivalent of 112.5 mmHg ($\pm 10\%$).

The pump will maintain the pre-set level. When operating, use the “+” and “-” pads (9) to increase and decrease the level. The factory preset level is 50.

4.3. Filling the chamber

Press the Fill Chamber pad to fill the pressure chamber initially, and throughout the surgery as needed. The chamber should never be over 1/3 full.

Note: numbers in (), refer to diagram on page 12.

4.4. Recommended initial settings.

Joint	Pressure level with Tourniquet	Pressure level without Tourniquet
Shoulder Joint		60
Acromioplasty		60
Knee Joint Hemarthrosis	30 50-60	65 50-60
Wrist	30	65
Elbow, Ankle	40	65
Hip		65

4.5. Pressure alarm

High pressure alarm.

If the high pressure safety level of 140 is exceeded, the pump stops and the red ALARM (6) LED lights up. The operator should lower the pressure to let the pump automatically restart. If the pump does not restart when the pressure drops below 140, restart the pump manually by turning the pump ON and OFF again with the POWER pad.

Low pressure alarm.

If the pressure falls below the low pressure safety level of 10, the red ALARM LED lights up, <Pres L> flashes on the LCD display, and the pump stops. Check that the luer-lock and the filter near the pressure transducer are in proper working order. If both are satisfactory press Fill Chamber to stop the alarm and Run/Stop pad to restart the unit.


Caution:

If water has accidentally entered the tubing leading to the pressure transducer, change the "One Day Set™" and restart the system.

Note: numbers in (), refer to diagram on page 12.

4.6. Pressure cycle mode.

This function increases the pressure by 60% for a two minute cycle. The pressure then returns automatically to its original setting. The pressure cycle mode is activated when bleeding occurs.

The pressure cycle mode is activated either by pressing the "Pressure" pad...

... or using the remote control, which is connected on the front panel of the pump...

... or by pressing on the red foot pedal Ref. 4171, which is connected on the rear panel of the pump (see on page 14, n° 16)

The remote control is not autoclavable and not soakable. A sterile disposable remote control cover must be used (see FMS Ref. 1087)

4.7. Maintenance: cleaning & decontamination

Cleaning the fms solo®:

At the end of the operating day, switch the fms solo® off, unplug the power cord and discard all tubing. Wipe off the instrument's sides, front panel and covers with a soft cloth that has been slightly moistened with a neutral - pH detergent. Wipe again using distilled water to rinse. **DO NOT SOAK.**

Do not use flammable agents when cleaning and disinfecting the fms solo®.

Foot pedal: cleaning and decontamination


After each operation, clean the foot pedal and the power cord with a neutral-pH. The foot pedal can be immersed, except the electrical connections which must remain dry.

DO NOT AUTOCLAVE.

Do not twist the cable when storing.

Do not unplug instrument by pulling on cable.

4.8. Interfaces cables

For installation of the cable or control unit, please refer to the instructions provided with the interface.

Shaver (Manufacturer)	Model / Class		fms hand control interface
FMS	TORNADO		Ref. 4120
LINVATEC	APEX (Class B)		Ref. 4107
LINVATEC	MICROCHOICE		Ref. 4103
LINVATEC	ADVANTAGE		Ref. 4104
DYONICS	EP-1		Ref. 4112
DYONICS	POWER		Ref. 4122
STRYKER	QUADRACUT		Ref. 4113
STRYKER	SE 5		Ref. 4127
STRYKER	TPS (ENDO PORT - Class B)		Ref. 4113
STRYKER	TPS 12K (TPS 1 or TPS 2 - Class B)		Ref. 4217*
STORZ	SL		Ref. 4138

When using the following shaver hand control interface, the activation of the pump functions is achieved with a 1 - way foot pedal. The activation of the shaver is achieved using the shaver manufacturer's hand control (button on shaver) or shaver foot pedal.

If your current shaver is not listed, contact fms.

*** Ref 4214 (110 V.) & 4217 (220 V.) have replaced Ref. 4114 & 4117. DO NOT USE Ref. 4114 & 4117. Call FMS for directions.**

Chapter 5. Troubleshooting

Problems	Cause	Possible solutions
The pump stops. Red ALARM LED lights up.	The HIGH PRESSURE safety has activated.	Lower the joint pressure. Press the Run/Stop pad to restart the pump.
	The LOW PRESSURE safety has activated: "Pres. L" blinks.	Verify that the luer-lock is tightly closed. Press the Fill Chamber to stop the alarm, and press the Run/Stop pad to restart the pump.
	The hydrophobic filter is wet.	Replace the "One Day Set™".
	There is no more irrigation fluid.	Spike the bag and press the Fill Chamber pad.
The roller pump heads do not rotate.	The power is on "O" or is on Stop.	Check that the main power switch on the back of the pump is on the "I" position, and that the Run/Stop pad green light is on.
There is no irrigation.	The stop-cock on the sheath or a clamp is closed, stopping outflow.	Check the inflow line from the bags to the joint. Open any closed clamps or stop-cocks.
The irrigation pump starts spinning frantically, emitting excessive noise.	There is no more irrigation fluid.	Replace the saline solution bags and/or check the clamps under the bags. Press Fill Chamber if chamber is lower than 1/3.

Problems	Cause	Possible solutions
The irrigation pump goes on and off erratically causing abnormal pressure fluctuations.	Water may have entered the pressure sensing line.	Change the “One Day Set™” tubing if water has entered the pressure sensing line.
There is too much water in the pressure chamber.	An air leak in the line may have developed.	Check that the luer-lock on the irrigation and pressure sensing tube are closed. If necessary, change the irrigation tubing and start again.
There is insufficient pressure in the system.	An irrigation problem exists.	Check the stop-cock on the arthroscope, and the clamps under the saline bags and under the pressure chamber.

Chapter 6. Product specification and ordering informations

Specifications

Composition

Ref: 4590

Dimensions

Height: 16 cm

Width: 39 cm

Depth: 30 cm

Weight: 11 kg

Performance specifications

Pressure levels : 20–140 Increments of 5

Precision $\pm 10\%$ at 0 flow rate

Electrical specifications

Input power: 120/230 V \sim

Frequency: 50–60 Hz

Current Draw: 500 VA

Fuses: 5 A Time Delay

Conforms to standard EN 60601-1

(Electromagnetic Compatibility)

Storage conditions

Temperatures: -10°C (50°F) to $+50^{\circ}\text{C}$ (104°F)

Operating conditions

Temperatures: $+10^{\circ}\text{C}$ to $+40^{\circ}\text{C}$

Relative humidity: 30 to 75%

Ordering Information

Equipment

	Ref
fms solo [®] (alone)	4590
1 - way foot pedal	4171
Remote control	8051

Disposables

“One Day Set [™] ” (ex. irrigation tubing)	4503
Patient set	1109
Irrigation tubing (ex. intermediary)	1102
fms Waste collection bag (10 liter)	1380
Remote control cover	1087

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fms solo[®] is a registered trademark of DePuy Mitek, Inc.
DePuy Mitek, Inc.
Future Medical Systems S.A.